

Infographic: Ultrasound QC for Base Unit, Transducer/Phantom, and HLD/Exam (Tropon) Traceability

QC-Track[®]

ATIRIX[®] Medical Systems

Always Inspection Ready[™]

Inspectors are looking more closely at Ultrasound QC than ever before, especially with US infection control a “TJC Top 10” concern. QC-Track helps you be *Always Inspection Ready*[™]. Here’s how.

Let’s start with the basics.

START

Ultrasound managers and infection control teams need to have three different types of tracking:

1. **Base Unit Device QC** - Can you prove that filters were regularly cleaned? Was the device ready for patient exams?
2. **Transducer/Phantom QC** - Can you prove that all transducers and base units were properly tested with a phantom?
3. **HLD/Exam Traceability** - Can you prove that transducers, especially intra-cavity or those in contact with non-intact skin, were properly cleaned prior to the patient’s exam? Will you be able to quickly identify and contact patients in the event of an outbreak?

...yet many facilities still track Ultrasound QC and traceability using spreadsheets and paper log books.

US teams have special tracking challenges, including **LOGISTICS**: large numbers of easily moved transducers, and **COMPLEXITY**: three different QC cycles for three different audiences.

QC-Track, with paperless QC workflows, makes it easy. How?

1

For Base Unit Device QC:

Audience:
Clinical Engineering, Physics Team

- QC-Track uses a scheduled workflow with checklist worksheets to track filters, overall unit cleanliness, monitor cleanliness, a visual check of transducers, and mechanical functions.
- Scheduled checklists can also be used to track the status of ancillary units, like a CleanShield[®] closet.

The result is consistent structure for US QC across your devices and locations.

2

For Transducer/Phantom QC:

Audience:
ACR, AIUM, Clinical Engineering, Physics Team

- QC-Track uses a scheduled workflow, with barcode tags on transducers and ultrasound base units, for transducer phantom QC.
- Reconciliation reports help ensure that all transducers are tested and give details on transducer and base unit testing history.

Note: Although annual per ACR CQC, clients have recently been advised by ACR inspectors to perform “Quality Control Tests” on at least a semi-annual basis, and some of our clients are doing this test monthly.

The result is high assurance that transducer phantom QC is being performed on all of your transducers.

3

For HLD/Exam Traceability:

Audience:
Hospital Infection Management, TJC

- Semi-critical devices, such as intra-cavity ultrasound transducers and those that come in contact with non-intact skin, are an area of concern for TJC. *See TJC Quick Safety Issue 33, displayed on back side.*
- QC-Track uses an on-demand workflow and QC Worksheets with adhesive barcode tags on transducers and cleaning units, such as Trophon, to track and link the patient exam with the prior cleaning record.
- The US team can track and report on the cleaning history of transducers and HLD cleaning systems, meeting the needs of both device QC tracking and TJC “Culture of Safety.”



The result is an easy, secure method for traceability and a powerful new tool for infection research.

THE RESULTS: An ultrasound QC process with all the benefits of paperless device QC

QC-Track provides the foundation for a comprehensive quality control program. Ultrasound managers, infection control teams, and inspectors have everything—reports, logs, and documents—at their fingertips. And you will appreciate the structure and consistency of paperless QC.

Ready? Contact Atirix now to learn more about QC-Track and enterprise imaging QC solutions for ultrasound and more. 877-273-1764 | atirix.com | QC-Answers@Atirix.com

Product or service names mentioned herein are the trademarks of their respective owners. Atirix products are covered by issued and pending patents, including US 8,478,610 “Medical Imaging Device Quality Control System and Method” and US 8,428,969 “System and Method for Tracking Imaging Quality” © 2024 Atirix Medical Systems, Inc. All Rights Reserved.

ATIRIX[®]
Medical Systems

Ultrasound Device QC Recommendations and Requirements

Here are samples of recommendations and standards from TJC, ACR, and IAC. They are only examples. Consult with your QC coordinator, biomedical engineer, medical physicist, and/or Infection Control Team to determine your facility's specific ultrasound device QC workflow requirements.

The Joint Commission - Top 10 Most Challenging Ambulatory Care Standards in 2019

09/02/2020

#1 - IC.02.02.01
 Reduce the risk of infections associated with medical equipment, devices and supplies.

Top 10 Most Challenging Ambulatory Care Standards in 2019 published 9 September 2020, accessed at www.jointcommission.org/resources/news-and-multimedia/blogs/ambulatory-buzz on 24 March 2021

IAC - Standards and Guidelines for Adult Echocardiography Accreditation

Section 2B: Adult Transesophageal Echocardiography Testing

STANDARD – Instrumentation

2.1B – Cardiac Ultrasound System

2.1.1B Ultrasound instruments utilized for transesophageal echocardiographic studies (TEEs) must include the echocardiographic imaging system requirements, as outlined in the [Section III: Adult Transesophageal Echocardiography Testing, STANDARD – Instrumentation](#).

2.2B – Transesophageal Ultrasound Transducer

2.2.1B Transesophageal ultrasound transducers must be those manufactured for the ultrasound system of the facility.

2.2.2B Transesophageal ultrasound transducers must incorporate multiphase imaging capabilities.

2.2.3B The manufacturer's guidelines must be followed for the appropriate care and cleaning of the TEE transducer and adhere to the appropriate infection control practices to prevent the transmission of disease. The structural and electrical integrity of the transducer must be documented in the routine TEE audit log. "Failed."

IAC Standards and Guidelines for Adult Echocardiography Accreditation published 27 April 2018, accessed at intersocietal.org/echo/seeking/echo_standards.htm on 24 March 2021

American College of Radiology - Ultrasound Accreditation Program Requirements

QC Test	Description
Physical and Mechanical Inspection	Assesses the mechanical integrity of the equipment, and the safety of patient and operator.
Image Uniformity and Artifact Survey	Identifies the presence of artifacts, often axial or lateral streaks in scans of uniform sections of a phantom. The use of "in-air" images (i.e., images acquired without the use of gel or phantom) may also be useful in detecting superficial artifacts (an in-air only approach to uniformity assessment for the routine QC test is reasonable and acceptable; a good implementation should provide baseline images of the reverberation/interference pattern for comparison to aid in detecting potential uniformity problems). All transducer ports on each scanner should be tested using at least 1 transducer. In the case of single probe, it is likely left plugged into the same port all the time, and other ports are not used. In this case, not testing the other ports would be acceptable for accreditation purposes.
Geometric Accuracy (Mechanically scanned transducers only)	Commonly involves use of the scanner calipers to measure known distances between test targets. Measurement is required only in the mechanically scanned directions.
Ultrasound Scanner Electronic Image Display Performance	Maintaining the performance of the image display is critical for providing the greatest diagnostic benefit of the scanner. They should also include worklist monitors only if used for primary interpretation (other than color analysis). Display characteristics that are evaluated may include gray scale response, presence of pixel defects, and overall image quality. These evaluations are typically performed using specialized test pattern images. See ACR/AAPM-SEM Technical Standard for Electronic Practice of Medical Imaging for additional information on tests and testing methods.
Primary Interpretation Display Performance	This test should be performed semiannually, or as judged appropriate based on the specific display technology, or prior QC testing data. Primary diagnostic displays may be electronic soft-copy displays on a PACS workstation or film. Display characteristics that are evaluated may include grayscale calibration, presence of pixel defects, and overall image quality. These evaluations are typically performed using specialized test pattern images. See ACR/AAPM-SEM Technical Standard for Electronic Practice of Medical Imaging for additional information on tests and testing methods.

Quality Control: Ultrasound and Breast Ultrasound (Revised 3-2-2021) published 2 March 2021, accessed at accreditation.support.acr.org/support/solutions/articles on 24 March 2021

The Joint Commission - Quick Safety Issue 33: Improperly sterilized or HLD equipment - a growing problem

05/22/2017

"Standard IC.02.02.01 requires organizations to reduce the risk of infections associated with medical equipment, devices and supplies. This standard is applicable to Joint Commission-accredited hospitals (HAP), critical access hospitals (CAH), ambulatory (AHC) and office-based surgery (OBS) facilities."

"According to reports to The Joint Commission's Office of Quality and Patient Safety, findings from non-complying organizations include:

- Lack of monitoring or documentation for sterilization or HLD of equipment, which makes it difficult to track the use of equipment on a specific patient, complicating the patient notification process when an outbreak occurs.
- Equipment is spread throughout the facility and may be processed or stored in numerous locations, making it difficult to track the equipment for documentation purposes."

Quick Safety 33: Improperly sterilized or HLD equipment—a growing problem, published 22 May 2017, accessed at www.jointcommission.org/resources/news-and-multimedia/newsletters/newsletters/quick-safety on 24 March 2021